



PRODUCT DATASHEET

Catalog No:	BSV-COV-PR-53
Product Name:	SARS-CoV-2 Nucleocapsid Protein (His Tag) (1 mg)
Description:	SARS-CoV-2 Nucleocapsid Protein with C-terminal polyhistidine tag.
Species:	2019-nCoV
2019-nCoV amino acids	2 – 419 (Gene Accession Number: YP_009724397)
Tag:	His-tag
Host:	<i>Escherichia coli</i>
Purity:	<p>>95% as determined by SDS-PAGE.</p> 

Predicted Molecular Mass:	<p>Calculated 46505.60 (experimental by EI at 46507.29 M+H+ ion)</p>  <p>Sample Name: Medical Comment: MNI-HT NCP Instrument Name: Instrument 1 Acquired Time: 17/07/2020 15:09:32 Data Filename: 19ms10.d ACQ Method: protein01_C18.m</p>
Homogeneity:	>95% determined by SDS-PAGE densitometry
Concentration:	5 mg/ml
Formulation:	PBS or HEPES buffered saline as required
Endotoxin:	< 1.0 ng/μg of protein (<1 EU/μg)
Storage and Shipping:	Store at -80°C. Avoid repeated freeze-thaw cycles. Shipped on dry ice.

FOR RESEARCH LABORATORY TEST USE ONLY!