

Nasopharyngeal Flocked Swab (BSV-905101)


Product description

- Swab is made by medical ABS plastic handle and flocking head.
- Application: Used for Nasopharyngeal sample collection, transportation and storage - selected for SARS-Cov-2.
- For medical use in Hospital, Centres for Disease Control, 3rd-party Testing Labs
- Shelf-life Lot: 1 year stored at room temperature
- Packaging NP swab: 30 pcs/box – 3000 pcs/carton - carton size: 74 x 37 x 24 – 3.8 kg



Product Specifications

Item Number	Total length	Break point	Flock Thickness	Tip length	Handle Diameter1	Handle Diameter2
905101	150 mm	100 mm	2 mm	22 mm	1.1 mm	2.5 mm



Product Certificates:

CE and DoC available

Ordering Information

Product Code	Description	Pack Size
BSV-905101-30	Nasopharyngeal Flocked Swabs (30 pcs)	30 pcs
BSV-905101-3000	Nasopharyngeal Flocked Swabs (3,000 pcs)	3,000 pcs