

PRODUCT DATA SHEET

HIGH BINDING CAPACITY

96 well plate

1. Description

Code	Plate format	Strips colour	Volume per well
MC01F-HB8	96 well solid plate	Clear	400 µl
MCW01F-HB8	96 well solid plate	White	400 µl
MCB01F-HB8	96 well solid plate	Black	400 µl

plate	Material	Polystyrene
	Well bottom	Flat
	Well shape	Round
Dimensions	SBS standard	

2. Features

Surface	High binding
Binding capacity	~ 400 ng IgG/cm ²
Uniformity	CV% < 5% with our standard test. For further details see Certificate of Quality
Packaging	Qty/ pk: 25
	Qty/ cs: 200
Storage	Room temperature
Shelf life	5 years
Other information	All lots are tested for uniformity and reproducibility Certificate of Quality is released for every lot